

Diaporthe perijuncta

Multiporator/Eppendorf Eporator®

Transformation Protocol

Protocol No. 4308 915.544 –09/2003

Microorganism	<i>Diaporthe perijuncta</i>
Cell type	Phytopathogenic fungus
Molecules injected	Virus RNA (<i>in vitro</i> -transcribed) from Diaporthe RNA virus
Growth medium	2 % potato dextrose agar (PDA)
Spheroplast solution	Chitinase-cellulase mixture, dissolved in 1 M MgSO ₄
Washing solution	1 M sorbitol; STC (1 M sorbitol, 50 mM Tris HCl (pH 8.0), 50 mM CaCl ₂)
Storage solution	Mixture of: STC, PTC, DMSO (80:20:1) [PTC: 40 % PEG 6000, 50 mM Tris HCl (pH 8), 50 mM CaCl ₂]
Electroporation solution	1 M sorbitol (ice-cold)
Outgrowth medium	1 g/l casein hydrolysate, 1 g/l yeast extract, 342 g/l sucrose, 16 g/l agar
Cuvette	1 mm gap width
Reference	Moleleki, N. et al • 2003 • Applied and Environmental Microbiology 69 • 3952-3956

Preparation of spheroplasts:

1. Suspend ten-day-old mycelium in a chitinase/cellulase solution overnight at room temperature.
2. Gravity filter the mycelium through a 120 µm pore size nylon mesh.
3. Add an equal volume of 1 M sorbitol (ice-cold) and centrifuge at 2,700 x g for 5 min at 4 °C.
4. Wash again with 1 M sorbitol, then resuspend pellet in 500 µl STC and centrifuge as above.
5. Resuspend spheroplasts in storage solution and use immediately or store at – 80 °C.

Electroporation of spheroplasts:

1. Resuspend spheroplasts in 85 µl of ice-cold 1 M sorbitol and place on ice.
2. Add RNase inhibitor (12 U) to a 15 µl transcription reaction mixture containing the *in vitro*-produced viral RNA transcripts.
3. Mix the viral transcripts with the spheroplast suspension and place on ice for 5 min.
4. Transfer mixture to a pre-chilled cuvette.
5. Wipe moisture from the cuvette and insert the cuvette into the device.
6. Electroporation:

Mode	Prokaryotes "O"
Voltage (V)	1,800 V
Time constant (τ)	5 ms
No. of pulses (n)	3 (manually)

7. Add 750 µl 1 M sorbitol (ice-cold) to the cuvette and place on ice for 10 min.
8. Pipette portions of 200 µl into 90-mm diameter petri dishes.
9. Add outgrowth medium at 48 °C to each petri dish and place plates on a laminar flow bench until the agar solidifies. Seal the petri dishes with parafilm and incubate overnight on a bench top at room temperature.
10. Incubate at 25 °C for 48 h, then for 5-10 days on a bench top at room temperature (20-25 °C).

eppendorf