



Cell type: hMSC

Performance of Eppendorf Cell Culture Consumables (TC treated)

Procedure: Human mesenchymal stem cells were seeded in appropriate culture medium (MSCGM™ BulletKit™). Cell Culture vessels were incubated at 37 °C with a humidified 5 % CO₂ atmosphere. Cell growth and morphology was monitored on day 4, 5 and 6 post-seeding.



Bars represent average count of two wells of an Eppendorf 6-Well Cell Culture Plate



Representative light microscope images (100X magnification)

The tissue culture (TC) treated surface of Eppendorf Cell Culture Consumables allows optimal attachment and homogeneous growth of **hMSC**.

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · 22331 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

MSCGM™ and BulletKit™ are trademarks of Lonza Group Ltd, Switzerland. Eppendorf® and the Eppendorf logo are registered trademarks of Eppendorf AG, Germany. U.S. Design Patents are listed on www.eppendorf.com/ip. All rights reserved, including graphics and images. Copyright © 2014 by Eppendorf AG.