

Eppendorf Certificate

Certificate of Quality – epT.I.P.S.[®] Motion Quality Assurance

All epT.I.P.S.[®] Motion pipette tips are single-use plastic tips which, in combination with all Eppendorf epMotion[®] automated liquid handling systems, form a perfectly functioning pipetting system.

Quality Assurance

All epT.I.P.S.[®] Motion pipette tips are manufactured in a controlled area and are subject to regular quality inspections within statistical process control.

The Eppendorf Quality Management System is certified according to the international standards ISO 9001, ISO 13485 and ISO 14001.

Composition

epT.I.P.S.[®] Motion pipette tips are made of virgin polypropylene of highest purity and quality. Our material suppliers do not use or intentionally incorporate the following agents into the materials Eppendorf uses for production:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyldodecylammonium salts (DiHEMDA))
- Plasticizers (softeners)

The used raw materials for epT.I.P.S.[®] Motion pipette tips comply with FDA regulations 21CFR § 177.1520 „Olefin Polymers“ and 21CFR § 178.2010 „Antioxidants and Stabilizers for Polymers“ for food containers. During production of epT.I.P.S.[®] Motion pipette tips no slip agents, biocides and plasticizers are used.

epT.I.P.S.[®] Motion pipette tips and all product components do not contain:

- Metallic dyes
- Latex
- Glucan
- Cellulose

Hamburg, January 2015



Dr. Carsten Buhlmann
Product Manager
Automation



Dr. Birgit Schreiber
Vice President
Quality Management & Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO
13485
Certified

ISO
14001
Certified