

Technical Data Sheet
Eppendorf Microplate VIS, 96/F-PS
English (EN)




Plate dimensions [mm]

Length			Width			Height	
l	l1	l2	w	w1	w2	h	fh
127.8	124.3	125.3	85.5	81.1	81.9	14.4	7.6

Well dimensions [mm]

Well dimensions [mm]			Positions of well center [mm]				
Diameter		Depth	Thickness	Position			
d	d1	h2	bt	P1	P2	P3	P4
7.0	6.8	11.3	1.4	14.4	9.0	11.2	9.0

Stack height [mm]

sh							
21.5							

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com