


Qualitätszertifikat/Certificate of Conformity

Quality certificate for rotor A-4-62 including four rectangular buckets

The rotor A-4-62 is an Eppendorf quality product which has been developed and produced in compliance with international standards. When used properly, the rotor has a service life of 7 years.

The bucket structures indicated in the pictures are technological surface structures which do not have any effect on the safe operation or the service life of the product. The safety of the rectangular bucket with this design has been verified by extensive calculations and endurance tests. The possibility of the bucket cracking during centrifugation due to these surface structures can therefore be excluded.



This does not apply in the case of subsequent damage which may be caused by improper use or exceeding the service life.

Proper use also includes cleaning, decontamination and maintenance of the rotor and the accessories. Further information can be found in the operating manual of the centrifuges 5810 and 5810 R.

Eppendorf accepts no warranty or liability for damage caused by incorrect or improper use.

Hamburg, November 2008

Eppendorf AG



Horst G. Schmidt
Director of Quality Management and Norms

eppendorf
In touch with life